

Quick Start Guide

PC Software and User Guide:

Download PC Software '**DigiCapture Pro**' and '**User Guide**' from www.download-toproview.com or scan below QR code. Also available from the attached CD.



App (TopScope Pro) :

iOS:

Search and install '**TopScope Pro**' from App Store.



iOS

Android :

Search and install '**TopScope Pro**' from play.google.com/store or download '**APK**' from www.download-topreview.com. Also available from the attached CD.



Android



APK(Scan browser)

How to connect the WiFi Microscope:

- 1: Turn on the WiFi setting on your iOS/Android Device, find '**WiFi_Microscope_xxx**' from the WiFi Device list.
- 2: Enter Password (the default password to connect the device is **123456789**)

Notice: More details find in the 'User Guide'.